


LMISC Install manual


25/10/2013

UMIN Center

1. Download compressed files.
2. Click the compressed files to start the self-extract process.
3. Select and open 「インストールモジュール」 folder.




4. Click setup program (setup.exe is the installer)




5. Click CRRedist2008_x86 to start installation.

6. Select and open 「上書きモジュール」 folder.



7. Select three files (LMISC.exe, lmisc_conf.xml, ODM.dll)



8. Copy (over-write) these three files into ¥LMISC folder.(¥LMISC folder is in the root directory of your HDD)